

CAMSCO SORBENT SELECTION CHART

| PHYSICAL PROPERTIES | | | | | SORBENT | | APPLICATIONS AND SELECTION GUIDE | | | | |
|----------------------------------|--------------------------------------|-------------------|------------------------|----------------------|---|-------------------------------------|---|---|---|---|--|
| Surface Area (m ² /g) | Packing Density (g/cm ³) | Maximum Temp (°C) | Conditioning Temp (°C) | Desorption Temp (°C) | (Equivalent or replaceable) | Mesh Sizes | Volatility Range (Carbon & BP) | Suitable Analytes | Features | Weakness/Caution | |
| 800 (1700) | N/A | 400 | 350 | 325 | Anasorb CMC (CMS) | Various | C3 ~ C4 -60 ~ 80°C | Developed for sampling very volatile organic compounds, freons, volatile halocarbons like methyl chloride and dichloromethane. Examples also include acetone, anaesthetic gases, propene and sulfurly fluoride (vicane) | High surface area, used for both thermal and solvent desorption. Slightly hydrophilic. Desorption efficiencies of polar compounds are higher than with charcoal when sampling VOCs. Comparable to Carbosieve SIII, similar to Carboxen 1000/1003 | Easily and irreversibly contaminated by high boilers – protect with front bed of weaker sorbent | |
| 980 | N/A | 350 | N/A Solvent | N/A Solvent | Anasorb 747 (Beaded active carbon) | 20/40 | C2 ~ C5 | A wide range of polar and nonpolar compounds, similar to charcoals but better for polar compounds. Examples include propene oxide, dichloromethane, methanol and a variety of ketones and acrylates. | High surface area, capacity for organic vapors similar to petroleum-based and coconut shell charcoal. Normally used with solvent desorption, but does not catalyze the breakdown of ketones on its surface | Easily and irreversibly contaminated by high boilers – protect with front bed of weaker sorbent | |
| 1060 | 0.45 | 400 | 350 | 330 | Carbosieve SII | 60/80 80/100 | C1 ~ C2 | Permanent gases (H ₂ , O ₂ , Ar, CO and CO ₂) and C1 ~ C2 hydrocarbons (methane, ethane, ethylene, acetylene) | High capacity / breakthrough volume for small molecules. Some hydrophilicity, low artifacts (<0.1ng) | Easily and irreversibly contaminated by high boilers – protect with front bed of weaker sorbent | |
| 820 | 0.61 | 400 | 350 | 330 | Carbosieve SIII | 60/80 | C2 ~ C4 -60 ~ 80°C | Permanent gases (H ₂ , O ₂ , Ar, CO and CO ₂) and C2 ~ C4 hydrocarbons, chloromethane | High capacity / breakthrough volume for small molecules. Moderately hydrophilic, low artifacts (<0.1ng). Comparable to Anasorb CMC, similar to Carboxen 1000/1003 | Low desorption efficiency for polar compounds. Less retentive capability than charcoal. Easily and irreversibly contaminated by high boilers. Retains more water than Carboxen 569 | |
| 1160 | 0.49 | 200 | 190 | 190 | Carbosieve G | Various | C1 ~ C3 | Permanent gases and C2-C3 hydrocarbons | High capacity / breakthrough volume for small molecules. Low back pressure | Easily and irreversibly contaminated by high boilers – protect with front bed of weaker sorbent | |
| 510 | 0.55 | 400 | 350 | 330 | Carboxen 563 (Ambersorb 340) | 20/45 | C3 ~ C5 50 ~ 200°C | C3 ~ C5 VOCs Similar to Carboxen 564 | Low back pressure. Highly hydrophobic – humidity proof. Preferred over Ambersorb XE-340 for higher capacity / breakthrough volume in VOC analysis | Low desorption efficiency for polar compounds. May produce sulfur compounds as artifacts, typically SO ₂ | |
| 400 | 0.59 | 400 | 350 | 330 | Carboxen 564 (Ambersorb 347) | 20/45 | C2 ~ C5 | C2 ~ C5 VOCs Similar to Carboxen 563, but less capacity for water | Low back pressure. Highly hydrophobic – humidity proof. Preferred over Ambersorb 347 for higher capacity / breakthrough volume in VOC analysis | Less sulfur content than Carboxen 563 | |
| 485 | 0.61 | 400 | 350 | 330 | Carboxen 569 | 20/45 | C2 ~ C5 | Similar to Carboxen 563 and 564, but higher capacity for organic molecules and less capacity for water | Closed micropores. Hydrophobic – humidity proof. No Ambersorb equivalent | | |
| 1200 | 0.52 | 400 | 350 | 300 | Carboxen 1000 (Carbosorb 572) (Purcosieve) | 40/60 60/80 80/100 (20/45) | C2 ~ C4 -60 ~ 80°C | Permanent gases (H ₂ , O ₂ , Ar, CO and CO ₂) and C2 ~ C4 hydrocarbons, i.e., vinyl chloride | Low artifacts but can produce sulfur compounds at high temp. Better desorption efficiency than S-III, significantly hydrophilic – do NOT use in humid conditions | Not as retentive as Carbosieve S-III. Easily and irreversibly contaminated by high boilers. May produce sulfur compounds as artifacts, typically SO ₂ . High artifacts (>10ng) | |
| 500 | 0.61 | 400 | 350 | 330 | Carboxen 1001 | 60/80 | C2 ~ C5 | Similar to Carboxen 569 | Similar to Carboxen 569 in strength and hydrophobicity | | |
| 1000 | 0.46 | 400 | 350 | 330 | Carboxen 1003 | 40/60 | C2 ~ C5 | Permanent gases (H ₂ , O ₂ , Ar, CO and CO ₂) and C2 ~ C5 hydrocarbons | Large surface area and hydrophobic. Similar to Carboxen 1000, Carbosieve S-III and Anasorb CMS, but retains even less water | Easily and irreversibly contaminated by high boilers – protect with front bed of weaker sorbent | |
| 1100 | 0.50 | 400 | 350 | 330 | Carboxen 1012 | | C4 ~ C6 | Larger micropore diameter than most Carboxen sorbents for slightly heavier analytes up to C6 | Inert. Also used for aqueous phase adsorption of organic compounds | Easily and irreversibly contaminated by high boilers – protect with front bed of weaker sorbent | |
| 75 | 0.40 | 400 | 350 | 330 | Carboxen 1016 | 60/80 | C3 ~ C9 | The weakest sorbent in the CMS family, extended range to C9 | Low surface area compared to other CMS sorbents – a very rare property in CMS family | Low breakthrough volume prohibits use as the final bed in a multi-sorbent design | |
| 675 (600) | 0.60 (0.62) | 400 | 350 | 330 | Carboxen 1018 (Carboxen 1021) | | C2 ~ C3 | Ethane, acetylene, acetaldehyde and other C2 ~ C3 | Breath sampling tubes exhaled breath samples. Inert and hydrophobic | | |
| 600 ~ 1000 | N/A | 400 | 350 | 300 | Molecular Sieve 5A | Various | -60 ~ 80°C | Permanent gases, and nitrous oxide | Significantly hydrophilic – do NOT use in humid conditions | Easily and irreversibly contaminated by high boilers, high artifacts (>10ng) | |
| 600 ~ 1000 | N/A | 350 | 330 | 300 | Molecular Sieve 13X | | -60 ~ 80°C | Permanent gases, and 1,3-butadiene | Significantly hydrophilic – do NOT use in humid conditions | Easily and irreversibly contaminated by high boilers, high artifacts (>10ng) | |
| 1200 | N/A | 400 | 375 | 350 | Unicarb (Spherocarb) | 60/80 | C3 ~ C8 -30 ~ 150°C | Both nonpolar and polar compounds. Especially very volatile, but specially large molecules, e.g. SF ₆ | Good batch-to-batch reproducibility, low artifacts (<0.1ng). Nonfriable and inert – suitable for labile compounds. Moderately hydrophilic | Easily contaminated by high boilers. Extremely expensive | |
| 560 | N/A | 400 | 350 | 350 | Carbograph 5TD | Various | C3 ~ C8 50 ~ 150°C | Light hydrocarbons C3 ~ C8 | High thermal stability. Low artifacts (<0.1ng). Hydrophobic – humidity proof | Some activity with labile compounds. Friable | |
| 100 (100) (100 ~ 200) | 0.35 (0.37) | 400 (400) (400) | 350 (350) | 325 (325) | Carbopack B (Carbotrap B) (Anasorb GCB1) (Carbograph 1) | 60/80 (20/40) (Various) | C5 ~ C12 >75°C (C5 ~ C12) (C5 ~ C12) | A wide range from medium to high volatility. Ketones, alcohols, and aldehydes (but not formaldehyde). Nonpolar within volatility range. Perfluorocarbon tracer gases. (GCB1 for both aliphatic and aromatic hydrocarbons) | High thermal stability. Low artifacts (<0.1ng). Low back pressure. Hydrophobic – humidity proof | Lower desorption efficiency than Tenax for higher molecular weight compounds when used for sampling in saturated atmospheres, i.e. during thermal stripping. Some activity with labile compounds. Friable | |
| 10 (10) (10 ~ 13) | 0.68 (0.68) | 400 | 350 | 325 | Carbopack C (Carbotrap C) (Anasorb GCB2) (Carbograph 2) | 60+ (20/40) (Various) (Various) | C8 ~ C20 | In general, Carbopack C extends the capability of Carbopack B to higher molecular weight analytes. Alkyl Benzenes and large aliphatics (C8 ~ C20). Heavy organics: PCBs (polychlorobiphenols), PNAs (Polynuclear aromatics) | High thermal stability. Low artifacts (<0.1ng). Low back pressure, low surface area. Hydrophobic – humidity proof | Some activity with labile compounds. Low surface area. Friable | |
| 5 (5) | 0.64 (0.69) | 400 | 350 | 325 | Carbopack F (Carbotrap F) | 60/80 (20/40) | C9 ~ C30 Esp. > C20 | In general, Carbopack F extends the capability of Carbopack B to higher molecular weight analytes. Lower surface area for trapping and efficiently releasing larger molecules, in the C9 ~ C30 range | High thermal stability. Very low surface area. Hydrophobic – humidity proof | Some activity with labile compounds. Very low surface area. Friable | |
| 240 (240) | 0.41 (0.43) | 400 | 400 | 400 | Carbopack X (Carbotrap X) | 40+ (60/80) (20/40) | C3 ~ C5 80 ~ 145°C | Light hydrocarbons: 1,3-Butadiene, benzene, toluene, and xylene | High thermal stability. Low artifacts (<0.1ng). Hydrophobic. Porous, bridges CMS's and GCB's in terms of sorbent strength | Some activity with labile compounds. Friable | |
| 24 (24) | 0.42 (0.45) | 400 | 350 | 325 | Carbopack Y (Carbotrap Y) | 60+ (20/40) | C12 ~ C20 | less volatile Hydrocarbons C12 ~ C20 | High thermal stability. Bridges Carbopack C and Carbopack B in terms of sorbent strength | Friable | |
| 220 | 0.18 | 400 | 350 | 325 | Carbopack Z | 60/80 (20/40) | C3 ~ C9 | light hydrocarbons C3 ~ C9 | High thermal stability | Friable | |
| <5 | 1.68 | 350 | 350 | 330 | Glass Beads | 60/80 | | Semi-volatiles, solids at room temperature | Thermally stable, inert, low surface area. Acts as a filter at tube inlet, segregating higher boiling compounds from more innocuous adsorbents | Suitable only for large molecules | |
| 750 | 0.76 | 200 | 190 | 180 | Silica Gel | 60/80 80/100 | Water, polar compounds | Low-boiling polar compounds, especially useful for separating chlorinated or sulfur compounds from matrices with hydrocarbon interferences | Within optimum range, good adsorption / desorption qualities. Either granular or beaded forms. Can be chemically modified to fit different applications | Retains water (adsorbs up to 40% of its own weight in moisture), little affinity for non-polar compounds | |
| 1070 (1200) | 0.57 | 400 | 190 | 180 | Coconut Charcoal (Anasorb CSC) | 20/40 60/80 | C2 ~ C5 -80 ~ 50°C | Wide variety of non-polar compounds | High capacity / breakthrough volume for low boilers. Greater retention capability (lower desorption efficiency) than Carbosieve SIII or Anasorb CMS | Higher capacity than Carbosieve SIII or Anasorb CMS. Tendency to retain water. May catalyze the breakdown of ketones | |
| 1050 | 0.50 | 400 | 190 | 180 | Petroleum (XC) Charcoal | 20/40 (20/40) | C2 ~ C5 -80 ~ 50°C | Charcoal derived from residue of petroleum products | Same as above | Same as above | |
| 330 | 1.02 | 190 | N/A Solvent | N/A Solvent | Amberlite XAD-2 (Supelapak-2) | 20/60 | Specific | Polyaromatic hydrocarbons, chlorinated pesticides, organoathio phosphates. Used to remove hydrophobic compounds up to 20,000 MW | Styrene/divinylbenzene copolymer. Moderate surface area; hydrophobic (dipole moment 0.3). Normally used with solvent desorption | Low temperature limit | |
| 725 | 1.02 | 150 | N/A Solvent | N/A Solvent | Amberlite XAD-4 | 20/60 | Specific | Used to remove small hydrophobic compounds and surfactants; widely used to remove chlorinated organics, organophosphorus pesticides, phenols, etc. | Styrene/divinylbenzene copolymer. Same as above, but larger surface area than XAD-2. Similarities to Anasorb 727 and Chromosorb 106 | Low temperature limit | |
| 700 ~ 800 | N/A | 250 | 250 | 225 | Anasorb 727 | 20/40 | C5 ~ C12 50 ~ 200°C | Developed for sampling reactive compounds thanks to its unreactiveness, but capable of a wide range of compounds similar to Chromosorb 106 | High surface area, highly inert. Extremely hydrophobic, all-purpose sorbent similar to XAD-4 and Chromosorb 106, thermal or solvent desorption | Low temperature limit. Special order item, limited mesh sizes available | |
| 350 | 0.29 | 250 | 250 | 225 | Chromosorb 102 | Various | C5 ~ C12 50 ~ 200°C | Wide range of VOCs, oxygenated compounds, halocarbons and chlorine-containing pesticides that has a BP >40°C (less volatile than methylene chloride). Some functional similarities to Amberlite XAD-2 | Styrene/Divinylbenzene (DVB) copolymer. Hydrophobic, inert. Some background at high temp (> 200°C). Note: The polarity of sorbents: Chromosorb 106 < Parapak Q < Chromosorb 102 < Parapak R < Chromosorb 105 < Parapak N < Chromosorb 101 < Parapak P < Chromosorb 103 < Chromosorb 104 | Low temperature limit. High artifacts (10ng) | |
| 700 ~ 800 | 0.28 | 250 | 250 | 225 | Chromosorb 106 | Various | C5 ~ C12 50 ~ 200°C | Low-boiling hydrocarbons, benzene, labile compounds, volatile oxygenated compounds. Functionally similar to Anasorb 727, more mesh sizes | Styrene/ DVB copolymer. Hydrophobic, inert. The least polar polymer in the Chromosorb family. Some functional similarities to Anasorb 727 and XAD-4 | Low temperature limit. High artifacts (10ng), batch-to-batch variations | |
| 400 ~ 500 | 0.30 | 250 | 225 | 225 | Chromosorb 107 | Various | Up to 150°C | Vinyl acetate, formaldehyde from water and acetylene from lower hydrocarbons. Sulfur compounds. Not recommended for glycols and amines | Crosslinked acrylic ester polymer. Hydrophobic and polar. Note: The polarity of Chromosorb 107/108 increases with their exposure to high temperature (~ 200°C, for example) | Low temperature limit | |
| 100 ~ 200 | 0.30 | 250 | 225 | 225 | Chromosorb 108 (Anasorb 708) | Various | Up to 150°C | Polar small molecules such as alcohols, water, aldehydes and glycols | Crosslinked acrylic. Hydrophobic and polar. See the note above | Low temperature limit | |
| 526 | 0.36 | 165 | 165 | 150 | HayeSep A | Various | C0 ~ C2 | Hydrogen, oxygen, argon, carbon oxides, nitric oxide, C2, hydrogen sulphide, and water | DVB polymer. Polarity 7/10, hydrophilic. Permanent CO gases at ambient temp; C2, hydrogen sulphide and water at higher temperature | Low temperature limit | |
| 608 | 0.33 | 190 | 190 | 180 | HayeSep B | Various | C1 ~ C2 | C1 and C2 amines, trace levels of ammonia and water | DVB and polyethyleneimine (PEI) copolymer. Polarity 8/10. Highly hydrophilic | Low temperature limit | |
| 442 (650) | 0.34 | 250 (250) | 250 | 225 | HayeSep C (Chromosorb 105) | Various | NOx, COx and SOx | Polar hydrocarbons such as hydrogen cyanide, ammonia, hydrogen sulphide and water | DVB and acrylonitrile (ACN) copolymer, Polarity 6/10, hydrophilic. Similar separation characteristics to Chromosorb 104 | Low temperature limit | |
| 795 | 0.33 | 290 | 190 | 180 | HayeSep D | Various | Up to 160°C | Low molecular weight compounds, esp. acetylene, halogen, and sulfur groups. CO and CO ₂ analysis | DVB polymer. Polarity 11/10, highly hydrophobic. Backup for Tenax where carbon based adsorbents are unsuitable. Retains low boiling compounds that breakthrough Tenax, esp. in saturated atmosphere | Low temperature limit. Moderate artifact level at upper temperature limit. High pressure drop | |
| 405 (300) | 0.36 (0.38) | 165 (190) | 165 (180) | 150 (160) | HayeSep N (Porapak N) | Various | C5 ~ C8 50 ~ 150°C | Volatile nitriles, e.g. acrylonitrile, acetonitrile, propionitrile. Pyridine, volatile alcohols, ethanol, methyl ethyl ketone | DVB and ethylene glycol dimethacrylate (EGDM) copolymer. Polarity 9/10, highly hydrophilic. Porapak polymers are in order of increasing polarity as types P, PS, Q, QS, R, S, N, and T | Low temperature limit. High artifacts ~ 10ng) | |
| 165 (150) | 0.42 (0.27) | 250 (250) | 250 (225) | 225 (225) | HayeSep P (Porapak P / PS) | Various | Up to 200°C | Esters, ethers, ketones, alcohols, Hydrocarbons, fatty acids, aldehydes and nitro compounds. Not recommended for amines and anilines | DVB/Styrene copolymer. Polarity 3/10. Porapak polymers are in order of increasing polarity as types P, PS, Q, QS, R, S, N, and T | Low temperature limit | |
| 582 (550) | 0.35 (0.34) | 275 (250) | 275 (225) | 250 (225) | HayeSep Q (Porapak Q / QS) | Various | C5 ~ C12 50 ~ 200°C | Some similarities to Chromosorb 106. Not for amines and anilines, not for nitric oxides | DVB polymer. Polarity 2/10, hydrophobic. Porapak polymers are in order of increasing polarity as types P, PS, Q, QS, R, S, N, and T | Low temperature limit. High artifacts, nitrated by nitrogen oxide gases | |
| 344 (525) | 0.32 (0.30) | 250 (250) | 250 (225) | 225 (225) | HayeSep R (Porapak R) | Various | Up to 200°C | Match the polarity of the analyte to the polarity of the sorbent, while covering the boiling temperature. Esters, ethers, nitriles and nitro compounds. Not recommended for glycols and amines | DVB/N-vinyl(2-pyrrolidone) (NV2P) copolymer. Polarity 5/10, both hydrophobic and hydrophilic. Porapak polymers are in order of increasing polarity as types P, PS, Q, QS, R, S, N, and T | Low temperature limit | |
| 583 (375) | 0.33 (0.35) | 250 (250) | 250 (225) | 225 (225) | HayeSep S (Porapak S) | Various | Up to 200°C | Match the polarity of the analyte to the polarity of the sorbent, while covering the boiling temperature (Amines, amides, alcohols, aldehydes, hydrazines and ketones. Not for acids, glycols and nitriles) | DVB/4-vinylpyridine (4VP) copolymer. Polarity 4/10, both hydrophobic and hydrophilic. Porapak polymers are in order of increasing polarity as types P, PS, Q, QS, R, S, N, and T | Low temperature limit. Reacts with nitroalkanes | |
| 250 (300) | 0.38 (0.43) | 165 (190) | 165 (180) | 150 (160) | HayeSep T (Porapak T) | Various | Up to 150°C | Match the polarity of the analyte to the polarity of the sorbent, while covering the boiling temperature | EGDM polymer. Polarity 10/10, highly hydrophilic. Porapak polymers are in order of increasing polarity as types P, PS, Q, QS, R, S, N, and T | Low temperature limit | |
| 35 | 0.28 | 350 | 320 | 300 | Tenax TA | 20/35 60/80 | C6 ~ C26 100 ~ 400°C | Aromatics (except benzene), nonpolar (BP>150°C), semi-volatile polar (BP>150°C). Note: Tenax TA has replaced Tenax GC for lower background signals | Poly (2,6-diphenyl-p-phenyleneoxide) polymer. Low surface area, within optimum range, will readily / efficiently release what it adsorbs and can be easily cleaned to a very low background. Inert and does not react with labile compounds. Low inherent artifacts (<1ng). Low affinity for water, hydrophobic | Low breakthrough volume. May form some artifacts when heated, reported sources are: CO ₂ , benzene, toluene, benzaldehyde, acetophenone | |
| 24 | 0.41 | 350 | 320 | 300 | Tenax GR | 20/35 60/80 | C7 ~ C30 100 ~ 450°C | Alkyl benzenes, PAH's (polycyclic aromatic hydrocarbons), PCB's (polychlorobiphenols) | Low surface area. Lower affinity for water than Tenax-TA. Chemical composition: 30% graphite carbon & 70% Tenax TA | Low breakthrough volume | |

CARBON MOLECULAR SIEVES (CMS)

SPECIALTY MATERIALS

GRAPHITIZED CARBON BLACK (GCB)

POROUS POLYMERS

